



Remember to clean & replace your equipment

Mask

Clean daily | Replace every 3 months

Headgear

Clean weekly | Replace every 6 months

Tubing

Clean weekly | Replace every 3 months

Disposable Fine Filters

Change every 2 weeks

Full Face Cushions

Replace every month

Nasal Cushions and Pillows

Replace every 2 weeks

Humidifier Chamber

Clean daily | Replace every 6 months

You must meet your insurance carrier's compliance requirements to order. Coverage will vary based on insurance plan. Copays and deductibles will apply.