

Medtronic-MiniMed 630G, 670G and 770G series Insulin Infusion Pumps

| | |
|--|---|
| Manufacturer name: | Medtronic |
| About the Recall: | The pumps were not pre-programmed with their basal rates or other verified settings (i.e., bolus wizard settings, sensor settings, etc.), which must be set up and saved on their pump prior to use. |
| Products(s) numbers/models that are affected: | MiniMed™ 630G: MMT-1715, MMT-1755, MMT-1754 MiniMed™ 670G: MMT-1780, MMT-1781, MMT-1782, MMT-1760, MMT-1761, MMT-1762, MMT-1740, MMT-1741, MMT-1742 MiniMed™ 770G: MMT-1880, MMT-1881, MMT-1882, MMT-1892, MMT-1891, MMT-1890 |
| URL Link(s) to Manufacturer Notice/Web Site: | https://www.medtronicdiabetes.com 1-800-646-4633, option 1 |
| Actions to be taken by customer/Patient: | Do not use your pump until you have consulted with your healthcare professional to determine the settings. |