

Manufacturer name:	nnd Medical Technology												
About the Recall:	This notice has been issued due to a manufacturing error resulting in a failure of the EasyOne FlowTube detection performed by the EasyOne Air. Consequently, measurements with the EasyOne Air cannot be performed using the EasyOne Flow Tubes of the affected lot												
Products(s) numbers/models that are affected:	<table border="1" data-bbox="467 772 1580 898"> <thead> <tr> <th data-bbox="467 772 602 821">MMS #</th> <th data-bbox="602 772 805 821">MFG Catalog #</th> <th data-bbox="805 772 1369 821">Description</th> <th data-bbox="1369 772 1580 821">Affected Lot(s)</th> </tr> </thead> <tbody> <tr> <td data-bbox="467 821 602 856">1055598</td> <td data-bbox="602 821 805 856">141-5050-50</td> <td data-bbox="805 821 1369 856">MOUTHPIECE, F/EASYONE AIR DISPP (50/CS)</td> <td data-bbox="1369 821 1580 856">AA240327</td> </tr> <tr> <td data-bbox="467 856 602 898">1055599</td> <td data-bbox="602 856 805 898">141-5050-200</td> <td data-bbox="805 856 1369 898">MOUTHPIECE, F/EASYONE AIR DISPP (200/CS)</td> <td data-bbox="1369 856 1580 898">AA240327</td> </tr> </tbody> </table>	MMS #	MFG Catalog #	Description	Affected Lot(s)	1055598	141-5050-50	MOUTHPIECE, F/EASYONE AIR DISPP (50/CS)	AA240327	1055599	141-5050-200	MOUTHPIECE, F/EASYONE AIR DISPP (200/CS)	AA240327
MMS #	MFG Catalog #	Description	Affected Lot(s)										
1055598	141-5050-50	MOUTHPIECE, F/EASYONE AIR DISPP (50/CS)	AA240327										
1055599	141-5050-200	MOUTHPIECE, F/EASYONE AIR DISPP (200/CS)	AA240327										
URL Link(s) to Manufacturer Notice/Web Site:	Nddmed.com												
Actions to be taken by customer/Patient:	Immediately discontinue use of any product matching the affected item(s) and lot number(s) listed above. If you have product affected by this notice, please contact nnd Medical Technologies for return, disposal and replacement of any product you may have.												