

<b>Manufacturer name:</b>	Medline Industries																				
<b>About the Recall:</b>	<p>Medline Industries, LP has issued a recall for specific item(s) and lot(s) of Maxorb Extra CMC/Alginate Wound Dressings. Medline’s contract manufacturer has initiated a recall due to defects on pouches that could compromise the device's ability to maintain a sterile barrier. Additionally, patches of burnt or cracked polyethylene have been identified on the inside face of the primary packaging that could also compromise sterility. A device with compromised sterility could result in an infection to the patient if the product is used. The contract manufacturer has indicated that if product has already been used on patients under a three-month time period, patients should be monitored for symptoms during routine follow up.</p>																				
<b>Products(s) numbers/models that are affected:</b>	<table border="1"> <thead> <tr> <th>MMS #</th> <th>MFG Catalog #</th> <th>Description</th> <th>Affected Lot(s)</th> </tr> </thead> <tbody> <tr> <td>907085</td> <td>MSC7048EP</td> <td>DRESSING, MAXORB X ALGINATE CAL 4X8 (5/BX 10BX/CS)</td> <td>W00070988</td> </tr> <tr> <td>476569</td> <td>MSC7044EP</td> <td>DRESSING, MAXORB XTRA CAL ALGINAT 4"X4" LF (10/BX 5BX/CS)</td> <td>W00070789</td> </tr> <tr> <td>738897</td> <td>MSC7048EP</td> <td>DRESSING, MAXORB X ALGINATE CAL 4X8 (50/BX 10BX/C</td> <td>W00070988</td> </tr> <tr> <td>721275</td> <td>MSC7044EP</td> <td>DRESSING, MAXORB EXT CAL/AGL 4X4 (50/CS)</td> <td>W00070789</td> </tr> </tbody> </table>	MMS #	MFG Catalog #	Description	Affected Lot(s)	907085	MSC7048EP	DRESSING, MAXORB X ALGINATE CAL 4X8 (5/BX 10BX/CS)	W00070988	476569	MSC7044EP	DRESSING, MAXORB XTRA CAL ALGINAT 4"X4" LF (10/BX 5BX/CS)	W00070789	738897	MSC7048EP	DRESSING, MAXORB X ALGINATE CAL 4X8 (50/BX 10BX/C	W00070988	721275	MSC7044EP	DRESSING, MAXORB EXT CAL/AGL 4X4 (50/CS)	W00070789
MMS #	MFG Catalog #	Description	Affected Lot(s)																		
907085	MSC7048EP	DRESSING, MAXORB X ALGINATE CAL 4X8 (5/BX 10BX/CS)	W00070988																		
476569	MSC7044EP	DRESSING, MAXORB XTRA CAL ALGINAT 4"X4" LF (10/BX 5BX/CS)	W00070789																		
738897	MSC7048EP	DRESSING, MAXORB X ALGINATE CAL 4X8 (50/BX 10BX/C	W00070988																		
721275	MSC7044EP	DRESSING, MAXORB EXT CAL/AGL 4X4 (50/CS)	W00070789																		
<b>URL Link(s) to Manufacturer Notice/Web Site:</b>	<p>Medline.com</p> <p>If you have any questions, contact the Recall Department at 866-359-1704 or recalls@medline.com.</p>																				
<b>Actions to be taken by customer/Patient:</b>	<p>Please check your supplies for any burn marks on the packaging, or any material found inside the sterile packaging. Should you find that you have product with burn marks or particles of material within the packaging discard the product and reach out to Advanced Medical Solutions for further directions.</p>																				