

<b>Manufacturer name:</b>	Hollister, Advanced medical Solutions								
<b>About the Recall:</b>	An Urgent Product Recall regarding one specific lot(s) of their CalciCare Calcium Alginate Dressing. This notice has been issued because AMS has become aware of a defect on primary packaging pouches, in which minor missing patches of Polyethylene have been detected. The defect on these pouches could compromise the device's ability to maintain a sterile barrier. In addition, patches of burnt or cracked polyethylene have been identified on the inside face of the primary packaging.								
<b>Products(s) numbers/models that are affected:</b>	<table border="1" data-bbox="386 911 1507 1003"> <thead> <tr> <th data-bbox="386 911 500 961">MMS #</th> <th data-bbox="506 911 685 961">MFG Catalog #</th> <th data-bbox="691 911 1318 961">Description</th> <th data-bbox="1325 911 1507 961">Affected Lot(s)</th> </tr> </thead> <tbody> <tr> <td data-bbox="386 970 500 1003">1012495</td> <td data-bbox="506 970 685 1003">529937R</td> <td data-bbox="691 970 1318 1003">DRESSING, ALCICARE CALCIUM ALGINATE 4X4.75 (10//BX)</td> <td data-bbox="1325 970 1507 1003">4D134</td> </tr> </tbody> </table>	MMS #	MFG Catalog #	Description	Affected Lot(s)	1012495	529937R	DRESSING, ALCICARE CALCIUM ALGINATE 4X4.75 (10//BX)	4D134
MMS #	MFG Catalog #	Description	Affected Lot(s)						
1012495	529937R	DRESSING, ALCICARE CALCIUM ALGINATE 4X4.75 (10//BX)	4D134						
<b>URL Link(s) to Manufacturer Notice/Web Site:</b>	<p data-bbox="360 1045 669 1087"><a href="http://www.hollister.com">www.hollister.com</a></p> <p data-bbox="360 1104 1562 1188">For questions regarding this notification, please contact Hollister at (888) 808-7456</p>								
<b>Actions to be taken by customer/Patient:</b>	<p data-bbox="360 1213 1578 1352">Immediately discontinue use of any product matching the affected item(s) and lot number(s) listed above. If you have no product matching the affected item(s) and lot number(s), no further action is needed.</p> <p data-bbox="360 1373 1578 1608">Please check your supplies for any burn marks on the packaging, or any material found inside the sterile packaging. Should you find that you have product with burn marks or particles of material within the packaging discard the product and reach out to Advanced Medical Solutions for further directions.</p>								