

<b>Manufacturer name:</b>	Advanced Medical Solutions																
<b>About the Recall:</b>	AdaptHealth has been made aware of a voluntary medical device recall from Advanced Medical Solutions regarding the following products; Silvercel Hydro Alginate, Silvercel Non-Adherent, Tegaderm Alginate, ActivHeal Alginate, ActivHeal Aquafibre, ActivHeal Non-Adhesive Foam, ActivHeal Non Adhesive Tracheostomy, ActivHeal PHMB Foam Non-Adhesive, ActivHeal Raponigel Ultra, Biatain Alginate, Hyalo 4 High Gelling Fire, Caldicare Reinforced Aquafibre, Nurocel Extra, Maxorb Extra. The defect on these pouches could compromise the device ability to maintain a sterile barrier. In addition, patches of burnt or cracked polyethylene have been identified on the inside face of the primary packaging.																
<b>Products(s) numbers/models that are affected:</b>	<table border="1" data-bbox="415 877 1523 1037"> <thead> <tr> <th data-bbox="415 877 521 932">MMS #</th> <th data-bbox="521 877 704 932">MFG Catalog #</th> <th data-bbox="704 877 1333 932">Description</th> <th data-bbox="1333 877 1523 932">Affected Lot(s)</th> </tr> </thead> <tbody> <tr> <td data-bbox="415 932 521 966">702428</td> <td data-bbox="521 932 704 966">800202</td> <td data-bbox="704 932 1333 966">DRESSING, SLVR ANTIMIC ALGISTR2X2 (10/BX) (SMS) SYSTGX</td> <td data-bbox="1333 932 1523 966">69223V002</td> </tr> <tr> <td data-bbox="415 966 521 999">746627</td> <td data-bbox="521 966 704 999">800202</td> <td data-bbox="704 966 1333 999">DRESSING, SLVRCEL ANTIMIC ALGISTR N/WVN 2X2(10/CT) SYSTGX</td> <td data-bbox="1333 966 1523 999">69223V002</td> </tr> <tr> <td data-bbox="415 999 521 1033">482617</td> <td data-bbox="521 999 704 1033">800202</td> <td data-bbox="704 999 1333 1033">DRESSING, SLVRCEL ANTIMIC ALGISTR N/WOVN 2X2 (50/C) SYSTGX</td> <td data-bbox="1333 999 1523 1033">69223V002</td> </tr> </tbody> </table>	MMS #	MFG Catalog #	Description	Affected Lot(s)	702428	800202	DRESSING, SLVR ANTIMIC ALGISTR2X2 (10/BX) (SMS) SYSTGX	69223V002	746627	800202	DRESSING, SLVRCEL ANTIMIC ALGISTR N/WVN 2X2(10/CT) SYSTGX	69223V002	482617	800202	DRESSING, SLVRCEL ANTIMIC ALGISTR N/WOVN 2X2 (50/C) SYSTGX	69223V002
MMS #	MFG Catalog #	Description	Affected Lot(s)														
702428	800202	DRESSING, SLVR ANTIMIC ALGISTR2X2 (10/BX) (SMS) SYSTGX	69223V002														
746627	800202	DRESSING, SLVRCEL ANTIMIC ALGISTR N/WVN 2X2(10/CT) SYSTGX	69223V002														
482617	800202	DRESSING, SLVRCEL ANTIMIC ALGISTR N/WOVN 2X2 (50/C) SYSTGX	69223V002														
<b>URL Link(s) to Manufacturer Notice/Web Site:</b>	admedsol.com																
<b>Actions to be taken by customer/Patient:</b>	Please check your supplies for any burn marks on the packaging, or any material found inside the sterile packaging. Should you find that you have product with burn marks or particles of material within the packaging discard the product and reach out to Advanced Medical Solutions for further directions.																