

URGENT: MEDICAL DEVICE RECALL LM5CA Rhythm Healthcare 5L Oxygen Concentrator

Manufacturer name:	Jiangsu Jumao X-Care Medical Equipment Co., Ltd.
About the Recall:	<p>AdaptHealth has been notified of a voluntary recall by Jiangsu Jumao X-Care Medical Equipment Co., Ltd. regarding the LM5CA Rhythm 5L Oxygen Concentrator.</p> <p>This is due to incidents of the device melting and fire during the use of the oxygen concentrators. The affected concentrators were manufactured on October 27, 2003.</p>
Products(s) numbers/models that are affected:	<p>Serial number range:</p> <p>23JBC003701-23JBC004440</p>
URL Link(s) to Manufacturer Notice/Web Site:	<p>https://www.jumaomedical.com/</p>
Actions to be taken by customer/Patient:	<p>Questions regarding this recall may be directed to Rhythm Healthcare directly via email at contactus@rhythmhc.com.</p>