

Manufacturer name:	AstraZeneca Pharmaceuticals LP																																													
About the Recall:	The AIRSUPRA (albuterol/budesonide) 90 mcg/80 msg Inhalation Aerosol canister for oral inhalation use. This voluntary recall is for specific sample lots listed below in which no more than 0.005% of the lot may contain units of AIRSUPRA in which the actuator mouthpiece is in the incorrect orientation resulting in non-delivery of dose.																																													
Products(s) numbers/models that are affected:	<table><tr><th>Lot</th><th>Expiry</th><th>Dates Distributed</th></tr><tr><td>6270044C00</td><td>10/2026</td><td>4/9/2024 - 11/14/2024</td></tr><tr><td>6270053C00</td><td>11/2026</td><td>9/12/2024 - 2/20/2025</td></tr><tr><td>6270045C00</td><td>11/2026</td><td>7/25/2024 – 3/3/2025</td></tr><tr><td>6270047C00</td><td>11/2026</td><td>8/12/2024 - 2/14/2025</td></tr><tr><td>6270056C00</td><td>11/2026</td><td>10/7/2024 - 2/28/2025</td></tr><tr><td>6270040D00</td><td>10/2026</td><td>7/1/2024 - 9/25/2024</td></tr><tr><td>6270064C00</td><td>12/2026</td><td>2/17/2025 - 4/25/2025</td></tr><tr><td>6270052C00</td><td>11/2026</td><td>8/8/2024 - 12/11/2024</td></tr><tr><td>6270034E00</td><td>10/2026</td><td>5/2/2024 - 10/31/2024</td></tr><tr><td>6270063E00</td><td>11/2026</td><td>1/27/2025 - 4/25/2025</td></tr><tr><td>6270071D00</td><td>12/2026</td><td>1/30/2025 - 4/17/2025</td></tr><tr><td>6270075D00</td><td>12/2026</td><td>3/24/2025 - 4/29/2025</td></tr><tr><td>6270107C00</td><td>9/2027</td><td>4/10/2025 - 4/29/2025</td></tr><tr><td>6270075F00</td><td>12/2026</td><td>3/31/2025 - 4/29/2025</td></tr></table>	Lot	Expiry	Dates Distributed	6270044C00	10/2026	4/9/2024 - 11/14/2024	6270053C00	11/2026	9/12/2024 - 2/20/2025	6270045C00	11/2026	7/25/2024 – 3/3/2025	6270047C00	11/2026	8/12/2024 - 2/14/2025	6270056C00	11/2026	10/7/2024 - 2/28/2025	6270040D00	10/2026	7/1/2024 - 9/25/2024	6270064C00	12/2026	2/17/2025 - 4/25/2025	6270052C00	11/2026	8/8/2024 - 12/11/2024	6270034E00	10/2026	5/2/2024 - 10/31/2024	6270063E00	11/2026	1/27/2025 - 4/25/2025	6270071D00	12/2026	1/30/2025 - 4/17/2025	6270075D00	12/2026	3/24/2025 - 4/29/2025	6270107C00	9/2027	4/10/2025 - 4/29/2025	6270075F00	12/2026	3/31/2025 - 4/29/2025
Lot	Expiry	Dates Distributed																																												
6270044C00	10/2026	4/9/2024 - 11/14/2024																																												
6270053C00	11/2026	9/12/2024 - 2/20/2025																																												
6270045C00	11/2026	7/25/2024 – 3/3/2025																																												
6270047C00	11/2026	8/12/2024 - 2/14/2025																																												
6270056C00	11/2026	10/7/2024 - 2/28/2025																																												
6270040D00	10/2026	7/1/2024 - 9/25/2024																																												
6270064C00	12/2026	2/17/2025 - 4/25/2025																																												
6270052C00	11/2026	8/8/2024 - 12/11/2024																																												
6270034E00	10/2026	5/2/2024 - 10/31/2024																																												
6270063E00	11/2026	1/27/2025 - 4/25/2025																																												
6270071D00	12/2026	1/30/2025 - 4/17/2025																																												
6270075D00	12/2026	3/24/2025 - 4/29/2025																																												
6270107C00	9/2027	4/10/2025 - 4/29/2025																																												
6270075F00	12/2026	3/31/2025 - 4/29/2025																																												
URL Link(s) to Manufacturer Notice/Web Site:	www.astrazeneca.com																																													
Actions to be taken by customer/Patient:	Please check your supply and should you have the affected lot number listed above please discontinue use, and reach out to Sedgwick Claims Management Services at 1-866-201-9064. Sedgwick Claims Management Services: email AS6579@sedgwick.com or phone 1-800-492-4319 for medical inquiries.																																													