

Manufacturer name:	Philips Respironics									
About the Recall:	When using certain circuits with Trilogy Evo, Trilogy Evo 02, Trilogy Evo Universal and Trilogy EV300 devices, Philips Respironics has become of a potential issue that can occur when using preset circuit calibration settings with certain Fisher and Paykel single limb circuits.									
Products(s) numbers/models that are affected:	<table><tr><th>F&P Circuit Part Number</th><th>Philips Part Number</th><th>Part Description</th></tr><tr><td>RT202</td><td>1141426</td><td>Single-limb inspiratory heated, 20 mm ID circuit</td></tr><tr><td>900MR810</td><td>1141549</td><td>Reusable single-limb Evatherm circuit</td></tr></table>	F&P Circuit Part Number	Philips Part Number	Part Description	RT202	1141426	Single-limb inspiratory heated, 20 mm ID circuit	900MR810	1141549	Reusable single-limb Evatherm circuit
F&P Circuit Part Number	Philips Part Number	Part Description								
RT202	1141426	Single-limb inspiratory heated, 20 mm ID circuit								
900MR810	1141549	Reusable single-limb Evatherm circuit								
URL Link(s) to Manufacturer Notice/Web Site:	www.usa.philips.com									
Actions to be taken by customer/Patient:	Prior to initialing mechanical ventilation, perform a circuit calibration anytime the circuit is connected. For further information or support, contact your local DME provider.									