

<b>Manufacturer name:</b>	ICU medical
<b>About the Recall:</b>	ICU Medical is issuing a medical device correction of a potential cassette bag leak on specific models and lots of 50 ml and 100 ml CADD Medication Cassette Reservoirs. Certain CADD 50 ml and 100 ml Medication Cassette Reservoirs may exhibit weakened medication bag welds along the bag perimeter inside the plastic housing. This could result in medication leakage along the bag perimeter.

**Products(s)  
numbers/models  
that are affected:**

21-7001-24	CADD™ Medication Cassette Reservoir with clamp and female Luer. Nonvented stopper included. 50mL	10610586022111	6022009	6037749	6022008		
21-7002-24	CADD™ Medication Cassette Reservoir with clamp and female Luer. Nonvented stopper included. 100mL	10610586022128	4468662	6005575	6063845	6037750	6037753
			6037752	6062688	6062691	6026828	6026831
			6037751	6062689	6026827	6026829	6022010
			6026830	6012451	6012433		
21-7300-24	CADD™ Yellow Medication Cassette Reservoir with flow stop, clamp and female Luer. Nonvented stopper included 100mL	10610586027192	6053989	6037756	6022016	6037757	6037758
			6037760	6062694	6037755	6062695	6026836
			6053986	6053985	6026835	6053987	6053990
			6062696	6053988	6037759	6022019	6012453
			6026832	6026833	6026834		
21-7301-24	CADD™ Medication Cassette Reservoir with flow stop, clamp and female Luer. Nonvented stopper included. 50mL	10610586027215	6005588	6022022	6085509	6005584	6022021
			6005583	6022023	6062697	6059296	6059295
			6059294	6022028	6022031	6022032	6026839
			6026837	6022030	6022026	6022024	6022025
21-7302-24	CADD™ Medication Cassette Reservoir with flow stop, clamp and female Luer. Nonvented stopper included. 100mL	10610586027239	6026861	6062705	6037781	6026872	6037778
			6053814	6054013	6037782	6026851	6037780
			6062714	6062702	6026868	6026863	6037774
			6062721	6054008	6026854	6022049	6022045
			6062718	6054010	6026867	6022040	6022047
			6054007	6056238	6026853	6022041	6022039
			6062723	6062712	6026875	6022053	6022034
			6062717	6026864	6022048	6022046	6012470
			6062716	6037775	6026869	6022051	6022033
			6062711	6037777	6026866	4468671	4468669
			6037779	6026870	6026873	6012472	6005591
			6062706	6037773	6026874	6012474	6012443
			6062701	6037772	6026862	6005593	6012477
			6037776	6005601	6026865	4468670	6026859
21-7600-24	NRFit CADD™ Yellow Medication Cassette Reservoir with NRFit™ connector with flow stop, yellow-striped tubing, clamp and female NRFit™ connector. Yellow stopper included 100mL	10610586044007	6054015				

**URL Link(s) to  
Manufacturer  
Notice/Web Site:**

[www.icumed.com](http://www.icumed.com)

**Actions to be  
taken by  
customer/Patient:**

Please review your supply of CADD Medication Cassette Reservoirs to confirm if any of the listed lots are within your possession. Should you have any supplies please contact ICU Medical for replacement.

ICU Medical may be reached at 1-800-258-5361 or via email [customerservice@icumed.com](mailto:customerservice@icumed.com)